

## ATP sulfurylase Yeast

**Description:** Adenosine 5" Triphosphate Sulfurylase Yeast Recombinant produced in E.Coli is a non-glycosylated, polypeptide chain containing 511 amino acids and having a Mw of 57.7 kDa. Adenosine 5" Triphosphate Sulfurylase Yeast Recombinant catalyzes the activation of sulfate by transferring sulfate to the adenine monophosphate moiety of ATP to form adenosine 5

Catalog #:ENPS-360

For research use only.

**Synonyms:** Sulfate adenyltransferase, EC 2.7.7.4, Sulfate adenylate transferase, SAT, ATP-sulfurylase, Methionine-requiring protein 3, ATPS.

**Source:** Escherichia Coli.

**Physical Appearance:** Sterile Filtered White lyophilized (freeze-dried) powder.

**Amino Acid Sequence:** MPAPHGGILQ DLIARDALKK NELLSEAQSS DILVWNLTPR  
QLCDIELILN GGFSPLTGFL NENDYSSVVT DSRLADGTLW TIPITLDVDE AFANQIKPDT  
RIALFQDDEI PIALTVQDV YKPNKTIEAE KVFRGDPEHP AISYLFNVAG DYYVGGSLA  
IQLPQHYDYP GLRKTPAQLR LEFQSRQWDR VVAFQTRNPM HRAHRELTVR AAREANAKVL  
IHPVVGLTKP GD

**Purity:** Greater than 95.0% as determined by SDS-PAGE.

### Formulation:

The ATP sulphurylase protein was lyophilized after dialysis against lyophilized from 10mM NaP buffer, 100mM NaCl, 10mM Lactose, 1% PEG pH 7.5 and 0.75mM DTT.

### Stability:

Lyophilized Sulfate adenylate transferase although stable at 4°C for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution ATPS should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

### Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### Solubility:

Spin vial before opening. Reconstitute ATP sulphurylase with 5mM NaP pH-7.5 & 0.75mM DTT at a concentration ranging from 0.1mg 1mg per ml. Can be diluted further into other aqueous buffers. pH range between 7.0 8.5 is best.

### Introduction:

ATP sulphurylase synthesizes adenosine 5-sulphatophosphate from ATP and inorganic SO<sub>4</sub><sup>2-</sup>. This is the first reaction of a two step sequence in the formation of active sulphate, adenosine 3-phosphate 5-sulphatophosphate, which is a sulphate donor for a wide variety of compounds and is also involved in the reduction of sulphate.

### Biological Activity:

12 Units/mg.

**To place an order, please [Click HERE](#).**