

CHIC2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 19kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CHIC2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BTL;

Catalog #: A7971

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26511

Isotype: IgG

Swiss Prot: Q9UKJ5

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the CHIC family of proteins. The encoded protein contains a cysteine-rich hydrophobic (CHIC) motif, and is localized to vesicular structures and the plasma membrane. This gene is associated with some cases of acute myeloid leukemia.

To place an order, please [Click HERE](#).