

MEX3C

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 69kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MEX3C

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RKHD2; BM-013; MEX-3C; RNF194;

Catalog #: A8393

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51320

Isotype: IgG

Swiss Prot: Q5U5Q3

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of a family of proteins with two K homology (KH) RNA-binding domains and a C-terminal RING-finger domain. The protein interacts with mRNA via the KH domains, and the protein shuttles between the nucleus and cytoplasm. Polymorphisms in this gene may contribute to hypertension.

To place an order, please [Click HERE](#).