

Phospho-RELB-S573

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 70kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S573 of human RELB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

kp

Synonym:

IREL; I-REL; REL-B;

To place an order, please [Click HERE](#).

Catalog #: AP0240

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5971

Isotype: IgG

Swiss Prot: Q01201

Purity: Affinity purification

For research use only.