

## MLF1

---

**Reactivity:** Human

**Tested applications:** WB IF

**Recommended Dilution:** WB 1:500 - 1:2000 IF 1:50 - 1:100

**Calculated MW:** 31kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human MLF1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:** A8012

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 4291

**Isotype:** IgG

**Swiss Prot:** P58340

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes an oncoprotein which is thought to play a role in the phenotypic determination of hemopoietic cells. Translocations between this gene and nucleophosmin have been associated with myelodysplastic syndrome and acute myeloid leukemia. Multiple transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*