

## CNTROB

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**101kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human CNTROB

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

LIP8; PP1221;

**Catalog #:**A8277

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**116840

**Isotype:**IgG

**Swiss Prot:**Q8N137

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a centrosomal protein that interacts with BRCA2, and is required for centriole duplication and cytokinesis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

*To place an order, please [Click HERE](#).*