

NID1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 136kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NID1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NID;

Catalog #: A8349

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4811

Isotype: IgG

Swiss Prot: P14543

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the nidogen family of basement membrane glycoproteins. The protein interacts with several other components of basement membranes, and may play a role in cell interactions with the extracellular matrix.

To place an order, please [Click HERE](#).