

MTIF3

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 32kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MTIF3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IF3mt;

Catalog #: A8524

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 219402

Isotype: IgG

Swiss Prot: Q9H2K0

Purity: Affinity purification

For research use only.

Background:

This gene encodes a translation initiation factor that is involved in mitochondrial protein synthesis. Polymorphism in this gene is associated with the onset of Parkinson's disease. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 5.

To place an order, please [Click HERE](#).