

Phospho-SYT2-T202

Reactivity: Human Mouse Rat

Tested applications: IF

Recommended Dilution: IF 1:100 - 1:200

Calculated MW: 46kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding T202 of human SYT2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Sytl;

Background:

Synaptotagmins, like SYT2, are integral membrane proteins of synaptic vesicles thought to serve as Ca(2+) sensors in the process of vesicular trafficking and exocytosis (Hilbush and Morgan, 1994

To place an order, please [Click HERE](#).

Catalog #: AP0460

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 127833

Isotype: IgG

Swiss Prot: Q8N9I0

Purity: Affinity purification

For research use only.