

SCF Human

Description: Stem Cell Factor Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 165 amino acids and having a molecular mass of 18409 Dalton. The SCF is purified by proprietary chromatographic techniques.

Synonyms: Kit ligand Precursor, C-kit ligand, SCF, Mast cell growth factor, MGF, SF, KL-1, Kitl, DKFZp686F2250.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Amino Acid Sequence: The sequence of the first five N-terminal amino acids was determined and was found to be Met-Glu-Gly-Ile-Cys.

Purity: Greater than 98.0% as determined by: (a) Analysis by SEC-HPLC. (b) Analysis by SDS-PAGE.

Formulation:

Lyophilized from a concentrated (1mg/ml) solution in water containing 0.02% NaHCO₃.

Stability:

Lyophilized KIT ligand although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution SCF should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

It is recommended to reconstitute the lyophilized Stem Cell Factor in sterile 18M-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

Introduction:

Stem cell factor / KIT ligand (SCF) is a cytokine which binds CD117(c-Kit). SCF is also known as "steel factor" or "c-kit ligand". SCF exists in two forms, cell surface bound SCF and soluble (or free) SCF. Soluble SCF is produced by the cleavage of surface bound SCF by metalloproteases. SCF is a growth factor important for the survival, proliferation, and differentiation of hematopoietic stem cells and other hematopoietic progenitor cells. One of its roles is to change the BFU-E (burst-forming unit-erythroid) cells, which are the earliest erythrocyte precursors in the erythrocytic series, into the CFU-E (colony-forming unit-erythroid).

Biological Activity:

The ED₅₀ as determined by the dose-dependant stimulation of Human TF-1 cells is < 2 ng/ml, corresponding to a Specific Activity of 500,000 IU/mg.

References:

1. Title: Direct binding of Cbl to Y568 and Y936 of the stem cell factor receptor/c-Kit is required for ligand-induced ubiquitination, internalization and degradation. Publication: Biochemical Journal

Immediate Publication. Published on 19 Jun 2006 as manuscript

BJ20060464. Link: http://www.biochemj.org/bj/imps_x/pdf/BJ20060464.pdf 2. Title: The D816V

Mutation of c-Kit Circumvents a Requirement for Src Family Kinases in c-Kit Signal

Transduction. Publication: March 5, 2009, doi: 10.1074/jbc.M808058200 April 24, 2009 The

Journal of Biological Chemistry, 284, 11039-11047.

Link: <http://www.jbc.org/content/284/17/11039.full> 3. Title: Structural Basis for c-KIT Inhibition by the

Suppressor of Cytokine Signaling 6 (SOCS6) Ubiquitin Ligase. Publication: First Published on

October 28, 2010, doi: 10.1074/jbc.M110.173526 January 7, 2011 The Journal of Biological

Chemistry, 286, 480-490. Link: <http://www.jbc.org/content/286/1/480.full>

Catalog #: CYPs-262

For research use only.

To place an order, please [Click HERE](#).