

PLAP Human

Description: Placental Alkaline Phosphatase Human Recombinant encoding 154-287 amino acids expressed in E.coli, shows a 41kDa band on SDS-PAGE (including GST tag). PLAP Human Recombinant is purified by proprietary chromatographic techniques.

Catalog #: ENPS-340

For research use only.

Synonyms: ALP, PLAP, Alkaline phosphatase placental type, EC 3.1.3.1, PLAP-1, Alkaline phosphatase Regan isozyme.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered liquid.

Formulation:

PLAP in 50mM Tris-Acetate, pH7.5, 1mM EDTA and 20% Glycerol.

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Placental alkaline phosphatase (PLAP) is a membrane-associated siaglycoprotein enzyme normally present at high concentration in syncytiotrophoblasts within the placenta during the third trimester of gestation. The expression of PLAP was originally thought to be restricted to term placenta but a human PLAP-like variant has been described which shares more than 85% homology with PLAP itself. PLAP is expressed only in normal term placenta, endocervix and fallopian tube and also in ovarian and proximal gastrointestinal tumors. It is also commonly expressed in germ cell tumors and more recently described in seminomas.

To place an order, please [Click HERE](#).