

## HEV ORF3

**Description:** The e.coli derived HEV protein is fused with beta-galactosidase at the N-Terminus and contains the HEV immunodominant ORF3 92-123 a.a.

**Catalog #:** HEPS-280

**Purity:** HEV ORF3 protein is >95% pure as determined by 10% PAGE (coomassie staining).

For research use only.

**Purification Method:**

HEV ORF3 protein was purified by proprietary chromatographic technique.

**Specificity:**

Immunoreactive with sera HEV-infected individuals.

**Formulation:**

20mM Tris-HCl, pH-8, 10mM B-ME and 8M urea.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Applications:**

HEV ORF3 antigen is suitable for ELISA and Western blots, excellent antigen for detection of HEV with minimal specificity problems.

**Introduction:**

Hepatitis E virus (HEV), the major etiologic agent of enterically transmitted non-A, non-B hepatitis worldwide, is a spherical, non-enveloped, single stranded RNA virus that is approximately 32 to 34 nm in diameter. HEV belongs to a genus of HEV-like viruses (unassigned genus). HEV has a single-stranded polyadenylated RNA genome of approximately 8 kb. Based on its physicochemical properties it is presumed to be a calici-like virus.

**Storage:**

HEV ORF3 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

**To place an order, please [Click HERE](#).**