

## HBsAg adw Mutant P-142-S

---

**Description:** Recombinant HBsAg adw antigen was mutated by replacing the Proline residue at position 142 with Serine.

**Catalog #:** HXPS-893

**Source:** Pichia pastoris.

For research use only.

**Physical Appearance:** Sterile Filtered clear solution.

**Purity:** Greater than 85.0% as determined by SDS-PAGE Coomassie staining.

**Formulation:**

Sterile Filtered solution containing 20mM Na<sub>2</sub>HPO<sub>4</sub>, 0.03M NaCl and 0.001% Thimerosal, pH-7.4.

**Stability:**

HBsAg Should be stored at 4°C.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Introduction:**

HBsAg is the surface antigen of the Hepatitis-B-Virus (HBV). The capsid of a virus has different surface proteins from the rest of the virus. The antigen is a protein that binds specifically on one of these surface proteins. It is commonly referred to as the Australian Antigen.

**To place an order, please [Click HERE](#).**