

YWHAQ Human

Description: YWHAQ Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 245 amino acids (1-245) and having a molecular mass of 27 kDa. YWHAQ is purified by proprietary chromatographic techniques.

Catalog #:PKPS-261

For research use only.

Synonyms: 14-3-3 theta, 14-3-3 tau, 14-3-3 T-cell, HS1, YWHAQ, 1C5, 14-3-3 Tau, Tyr-3/Trp- 5 Monooxygenase Activation Protein Theta.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MEKTELIQKA KLA EQAERYD DMATCMKAVT EQGAELSNEE
RNLSSVAYKN VVGRRSAWRVISSIEQKTD TSDKKLQLIK DYREKVESEL RSICTTVLEL
LDKYLIANAT NPESKVFYLMKGDYFRYLA EVACGDDRKQ TIDNSQGAYQ EAFDISKKEM
QPTHPIRLGL ALNFSVFYIEILNNPELACT LAKTAFDEAI AELDTLNEDS YKDSTLIMQL
LRDNLTLWTS DSAGE

Purity: Greater than 95.0% as determined by (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.

Formulation:

YWHAQ solution containing 20mM Tris 7.5.

Stability:

YWHAQ Human Recombinant although stable at 4°C for 1 week, should be stored desiccated below -18°C. Please prevent freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

The 14-3-3 family of proteins plays a key regulatory role in signal transduction, checkpoint control, apoptotic and nutrient-sensing pathways. 14-3-3 proteins are highly conserved and ubiquitously expressed. There are at least seven isoforms, γ , δ , ϵ , ζ , η , θ , and τ that have been identified in mammals. The 14-3-3 tau, a subtype of the 14-3-3 family of proteins, was found in T Cells, brain and testes. This 14-3-3 tau is upregulated in patients with amyotrophic lateral sclerosis.

To place an order, please [Click HERE](#).