

## AFP Human

---

**Description:** Human alpha-fetoprotein purified from human cord serum.

**Catalog #:** PRPS-413

**Synonyms:** Alpha-fetoprotein, Alpha-fetoglobulin, Alpha-1-fetoprotein, AFP, FETA, HPAFP.

For research use only.

**Source:** Human cord serum.

**Physical Appearance:** Sterile Filtered clear solution.

**Formulation:**

Alpha Fetoprotein protein at 2.5mg/ml 0.1M PBS buffer pH 7.4 containing 0.1% NaN<sub>3</sub>.

**Stability:**

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Introduction:**

AFP is normally synthesized in the liver, intestinal tract, and yolk sac of the fetus. Antibody to AFP has been shown to be useful in detecting hepatocellular carcinomas (HCC) and germ cell neoplasms, especially yolk sac tumors.

**To place an order, please [Click HERE](#).**